



Certificates

FDA Registered Company

FDA 510(k) Registered Products

ISO & MDSAP Certified

We Private Label



Products Available In Bulk Packing By Inquiry

**FDA
510(k)**



Insulin Syringes	1
U- 40 Insulin Syringes	2
Tuberculin Syringes	3
1mL Syringes	3
3mL Syringes	4
5mL Syringes	5
10mL Syringes	5
12mL Curved Tip Syringes	5
20mL Syringes	6
30mL Syringes	6
50mL Syringes	6
Aluminum Hub Needles : Hard Pack	7
Aluminum Hub Blunt Needles	7
Hypodermic Needles	8
Safety Hypodermic Needles	9
1mL Syringe with Safety Needle	10
3mL Syringe with Safety Needle	10
Vaculet Blood Collection Sets	11
Scalp Vein (Butterfly) Sets	11
Multi-Sample Needles	12
Multi-Sample Luer Adapters	12
Multi-Sample Luer Holders	12
Multi-Sample Luer Lock Holders with Pre-Attached	12
Injection Plugs	12
Spinal Needles	13
IV Catheters	14
Surgical Blades	15
Disposable Scalpels	15
IV Administration Sets	16
IV Extension Sets	17
Bulk Non-Sterile Syringes	18
Product Specifications	19
Notes	20

INSULIN SYRINGES

Comfort point[®]

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Permanently attached needle
- ✓ Reduced dead space for reduced Insulin waste
- ✓ Smooth plunger action for greater filling accuracy
- ✓ Auto ground needle tip for greater comfort



0.3mL LO-DOSE INSULIN SYRINGE

with permanently attached needle

Catalog No.	Description
26018	U-100 29G x 1/2"
26014	U-100 30G x 5/16"
26010	U-100 31G x 5/16"
10/bag 100/box 500/case	

1/2mL LO-DOSE INSULIN SYRINGE

with permanently attached needle

Catalog No.	Description
26026	U-100 28G x 1/2"
26028	U-100 29G x 1/2"
26015	U-100 30G x 5/16"
26011	U-100 31G x 5/16"
10/bag 100/box 500/case	

1mL INSULIN SYRINGE

with permanently attached needle


Catalog No.	Description
26027	U-100 28G x 1/2"
26029	U-100 29G x 1/2"
26016	U-100 30G x 5/16"
26012	U-100 31G x 5/16"
10/bag 100/box 500/case	

Comfort point[®]

U - 40 INSULIN SYRINGES



- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Permanently attached needle
- ✓ Reduced dead space for reduced Insulin waste
- ✓ Smooth plunger action for greater filling accuracy
- ✓ Auto ground needle tip for greater comfort



U-40 INSULIN SYRINGE

with permanently attached needle

Catalog No.	Description
26036	U-40 1/2mL 29G x 1/2"
26037	U-40 1mL 29G x 1/2"

10/bag 100/box 500/case



EXEL U-40 *Comfort point*
Sterile Insulin Syringes

Comfort point a name recognized for quality at affordable pricing.
Ask your EXEL distributor for EXEL's U-40 Comfort Point Insulin Syringes ...
If it's not Comfort Point U-40 it's not EXEL. ... and start saving now.
EXEL is your proven partner in animal healthcare.

TUBERCULIN SYRINGES

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Smooth plunger action for greater filling accuracy
- ✓ Auto ground needle tip for greater comfort
- ✓ Reduced dead space for reduced medication waste



1mL TUBERCULIN SYRINGE

with detachable needle, Luer-slip

Catalog No.	Description
26040*	1mL 27G x 1/2"
26042*	1mL 25G x 1"
26043*	1mL 26G x 1/2"
26044*	1mL 25G x 5/8"
26045*	1mL 26G x 3/8"

100/box 1,000/case

1mL ALLERGY SYRINGE

with permanently attached needle

Catalog No.	Description
26061	1mL 27G x 1/2"
26063	1mL 27G x 3/8"

25/Tray 1000pc/Case



1mL TUBERCULIN SYRINGE ONLY

Catalog No.	Description
26048*	1mL Slip Tip with cap
26049*	1mL Lock Tip with cap
26050	1mL Lock Tip

100/box 1,000/case



Luer Slip

Luer Lock



*Low Dead Space Plunger

1mL ZERO DEAD SPACE TUBERCULIN SYRINGE

with permanently attached needle

Catalog No.	Description
26046	1mL 25G x 5/8"

10/bag 100/box 500/case

SINGLE UNIT BLISTERS

Catalog No.	Description
26051	1cc 21Gx1"
26066	1cc 25Gx1"
26067	1cc 23Gx1"

100/box 1,000/case



3mL SYRINGES

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Positive plunger stop to prevent spillage when drawing
- ✓ Smooth plunger motion
- ✓ Polypropylene construction



LUER-LOCK TIP

3mL Syringe / Needle Combination

Catalog No.	Description	Hub Color	
26100	25G x 5/8"	Orange	■
26101	23G x 1"	Deep Blue	■
26102	22G x 1"	Black	■
26103	22G x 1 1/4"	Black	■
26104	22G x 1 1/2"	Black	■
26105	21G x 1"	Deep Green	■
26107	21G x 1 1/2"	Deep Green	■
26108	20G x 1"	Yellow	■
26109	20G x 1 1/2"	Yellow	■
26110	18G x 1 1/2"	Pink	■
26111	25G x 1"	Orange	■
26112	25G x 1 1/2"	Orange	■
26115	22G x 3/4"	Black	■
26117	23G x 1 1/2"	Deep Blue	■
26129	27G x 1 1/4"	Medium Grey	■
100/box 1,000/case			

LUER-SLIP TIP

3mL Syringe / Needle Combination

Catalog No.	Description	Hub Color	
26601	22G x 3/4"	Black	■
26602	25G x 5/8"	Orange	■
26603	23G x 1"	Deep Blue	■
26605	22G x 1"	Black	■
26608	21G x 1"	Deep Green	■
26611	20G x 1"	Yellow	■

100/box 1,000/case

3mL SYRINGE only

Catalog No.	Description
26200	Luer Lock with cap
26201	Luer Slip with cap
100/box 1,000/case	

5mL SYRINGES

10mL SYRINGES

CURVED TIP SYRINGES

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ With Cap
- ✓ Easy to read large bold print
- ✓ Positive plunger stop to prevent spillage when drawing
- ✓ Smooth plunger motion
- ✓ Crystal Clear Barrels
- ✓ 5mL syringes are graduated up to 6mL of total graduations
- ✓ 10mL syringes are graduated up to 12mL of total graduations



26210



26213



26230



26251



26267

5mL SYRINGE/NEEDLE COMBINATION LUER-LOCK TIP

Catalog No.	Description	Hub Color
26210	22G x 1"	Black
26211	22G x 1½"	Black
26212	21G x 1"	Deep Green
26213	21G x 1½"	Deep Green
26214	20G x 1"	Yellow
26218	21G x 1¼"	Deep Green
26219	22G x 1¼"	Black

100/box 800/case

5mL SYRINGE only

Catalog No.	Description
26230	Luer Lock with cap
26231	Luer Slip with cap

100/box 800/case

10mL SYRINGE/NEEDLE COMBINATION LUER-LOCK TIP

Catalog No.	Description	Hub Color
26250	22G x 1"	Black
26251	22G x 1 ½"	Black
26252	21G x 1"	Deep Green
26253	21G x 1 ½"	Deep Green
26254	20G x 1"	Yellow
26255	20G x 1 ½"	Yellow
26257	18G x 1"	Pink

100/box 800/case

10mL SYRINGE only

Catalog No.	Description
26265	Luer Lock with cap
26266	Luer Slip with cap

100/box 800/case

12mL CURVED TIP SYRINGE

Catalog No.	Description
26267	Curved Tip

50/box 500/case (non-sterile)


20mL SYRINGES

30mL SYRINGES

50mL SYRINGES

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Easy to read large bold print
- ✓ Positive plunger stop to prevent spillage when drawing

- ✓ Smooth plunger motion
- ✓ Crystal Clear barrels
- ✓ Polypropylene construction
- ✓ 20mL syringes are graduated up to 25mL of total graduations
- ✓ 30mL syringes are graduated up to 35mL of total graduations
- ✓ 50mL syringes are graduated up to 60mL of total graduations



20mL SYRINGE only

Catalog No.	Description
26280	Luer Lock with cap
26281	Luer Slip Ecc with cap
50/box 500/case	

30mL SYRINGE only

Catalog No.	Description
26290	Luer Lock with cap
26291	Luer Slip Ecc with cap
26292	Cath Tip Ecc
50/box 250/case	

50mL SYRINGE only

Catalog No.	Description
26300	Luer Lock with cap
26301	Luer Slip Ecc with cap
26304	Catheter Tip with cap
25/box 150/case	

ALUMINUM HUB NEEDLES

ALUMINUM HUB BLUNT NEEDLES

- ✓ Designed for veterinary use
- ✓ Color - coded for easy identification
- ✓ Sterile, Non-Toxic, Non-Pyrogenic
- ✓ Each needle is individually packed and sealed with a tamper evident heat stake
- ✓ Not for Human Use

ALUMINUM HUB HYPODERMIC NEEDLES - HARD PACK

Catalog No.	Description	Hub Color	
26450	14G x 1"	Tan	■
26452	14G x 1 ½"	Tan	■
26454	14G x 2"	Tan	■
26455B	16G x ⅝"	Gray	■
26456	16G x ¾"	Gray	■
26457	16G x 1"	Gray	■
26458	16G x 1 ½"	Gray	■
26461	18G x ¾"	Green	■
26462	18G x 1"	Green	■
26463B	18G x ⅝"	Green	■
26464	18G x 1 ½"	Green	■
26467	20G x 1"	Pink	■
100/box 1000/case			



ALUMINUM HUB BLUNT NEEDLES

Catalog No.	Description	Hub Colors	
26471	15G x 1 ½"	Clear	□
26472	17G x 1 ½"	Pale Green	■
26473	18G x 1"	Green	■
26475	20G x 1 ½"	Pink	■
25/box, 100/case			



HYPODERMIC NEEDLES



- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ SUS - 304 stainless steel
- ✓ Needles and boxes are color-coded for easy identification
- ✓ Each needle is individually packaged, sterilized and ready for use

SPECIALTY USE HYPODERMIC NEEDLES

Catalog No.	Description	Hub Color	
26437	30G x ½"	Yellow	■
26439	30G x 1"	Yellow	■
100/box 2,000/case			
26448	28G x ¾"	Green Blue	■
100/box 1,000/case			
26391	31G x ½"	Orange	■
26392	32G x ½"	Green	■
26392S	32G x ½"	Green(soft pack)	■
26393	33G x ½"	Clear	□
100/box 1,000/case			

NEEDLE LENGTH



Inches	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"
Millimeters	13mm	16mm	19mm	25mm	32mm	38mm	50mm

HYPODERMIC NEEDLE REGULAR BEVEL


Description	Catalog No.	Hub Color	
26400	27G x ½"	Medium Grey	■
26402	26G x ½"	Brown	■
26403	25G x 5/8"	Orange	■
26405	25G x 1"	Orange	■
26406	25G x 1 ½"	Orange	■
26407	23G x ¾"	Blue	■
26408	23G x 1"	Blue	■
26409	23G x 1 ½"	Blue	■
26410	22G x ¾"	Black	■
26411	22G x 1"	Black	■
26412	22G x 1 ½"	Black	■
26414	21G x 1"	Green	■
26417	20G x 1"	Yellow	■
26418	20G x 1 ½"	Yellow	■
26419	18G x 1"	Pink	■
26420	18G x 1 ½"	Pink	■
26424	19G x 1"	Cream	■
26426	27G x 1 ½"	Medium Grey	■
26427	27G x 1 ¼"	Medium Grey	■
26430	20G x ¾"	Yellow	■
26431	20G x ½"	Yellow	■
26438	19G x 1 ½"	Cream	■
100/box 2,000/case			

SecureTouch[®]

SAFETY HYPODERMIC NEEDLES



OR



- ✓ This device does not contain natural rubber latex.
- ✓ Sterile, Non-Toxic, Non-Pyrogenic.
- ✓ One-Handed Activation
- ✓ No Solid surface required
- ✓ Simple finger stroke activation
- ✓ Protective arm keeps fingers behind the needle
- ✓ An audible “click” as well as tactile feel and a visual indicator confirm safety mechanism engagement
- ✓ Needle remains locked inside activated cover
- ✓ Fits any brand syringe.
- ✓ Safety shield can be pushed back comparative to the needle, allowing a clear view of the injection site.

SAFETY HYPODERMIC NEEDLES

Catalog No.	Description	Hub Color	
27403	25G x 5/8"	Orange	■
27405	25G x 1"	Orange	■
27406	25G x 1 1/2"	Orange	■
27408	23G x 1"	Deep Blue	■
27417	20G x 1"	Yellow	■
27420	18G x 1 1/2"	Pink	■
27426	27G x 1 1/2"	Medium Grey	■
27427	27G x 1 1/4"	Medium Grey	■
27437	30G x 1/2"	Yellow	■

100/box 1000/case

* Made to order; lead times may apply.

SYRINGE WITH SAFETY HYPODERMIC NEEDLES

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile, Non-Toxic, Non-Pyrogenic.
- ✓ One-Handed Activation
- ✓ No Solid surface required
- ✓ Simple finger stroke activation
- ✓ Protective arm keeps fingers behind the needle
- ✓ An audible “click” as well as tactile feel and a visual indicator confirm safety mechanism engagement
- ✓ Needle remains locked inside activated cover
- ✓ Safety shield can be pushed back comparative to the needle, allowing a clear view of the injection site.

SecureTouch[®]

SYRINGE WITH SAFETY NEEDLE




3mL LL SYRINGE


Catalog No.	Description	Hub Color	Catalog No.	Description	Hub Color
27100	25G x 5/8"	Orange	27107	21G x 1 1/2"	Deep Green
27101	23G x 1"	Deep Blue	27108	20G x 1"	Yellow
27102	22G x 1"	Black	27110	18G x 1 1/2"	Pink
27104	22G x 1 1/2"	Black	27111	25G x 1"	Orange
27105	21G x 1"	Deep Green	27112	25G x 1 1/2"	Orange

50/box 400/case

- 3mL max dead space without needle 0.029mL (29µL)



OR



1mL LL SYRINGE

Catalog No. Description Hub Color

27042	25G X 1"	Orange
27044	25G X 5/8"	Orange
27047	23G X 1"	Deep Blue

50/box 300/case


- 1mL max dead space without needle 0.019mL (19µL)

* Made to order; lead times may apply.

BLOOD COLLECTION SETS

SCALP VEIN SETS

- ✓ This product does not contain Natural Rubber Latex.
- ✓ Individually packed.
- ✓ Thin wall stainless steel cannula, siliconized for easy insertion.
- ✓ Non-Pyrogenic
- ✓ Flexible wings hold the needle in position and help to prevent inadvertent pull out.
- ✓ Color-coded for easier identification.



SCALP VEIN SET PSV / BUTTERFLY / SMALL VEIN INFUSION

Catalog No.	Description	Tube	Hub Color	
26702	19G x 3/4"	12"	Cream	
26703	20G x 3/4"	12"	Yellow	
26704	21G x 3/4"	12"	Green	
26706	23G x 3/4"	12"	Deep Blue	
26708	25G x 3/4"	12"	Orange	
26709	27G x 3/4"	12"	Medium Grey	
50/box 500/case				

BLOOD COLLECTION SETS

Catalog No.	Description	Tube	Hub Color	
26764	21G x 3/4"	12"	Green	
26766	23G x 3/4"	12"	Deep Blue	
26768	25G x 3/4"	12"	Orange	
50/box 200/case				



SecureTouch®

SAFETY SCALP VEIN SET SAFETY PSV / BUTTERFLY / SMALL VEIN INFUSION

Catalog No.	Description	Tube	Hub Color	
27704	21G x 3/4"	12"	Green	
27706	23G x 3/4"	12"	Deep Blue	
27708	25G x 3/4"	12"	Orange	
*NON DEHP 50/box 500/case				



MULTI-SAMPLE NEEDLES

MULTI-SAMPLE LUER ADAPTERS

MULTI-SAMPLE LUER HOLDERS

INJECTION PLUG

- ✓ This product does not contain Natural Rubber Latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Multi-sample needle permit several samples to be taken with a single puncture
- ✓ Sharp & smooth edges make penetration painless for patients and provides easy connection to rubber stoppers

MULTI-SAMPLE (BLOOD DRAW) NEEDLES

Catalog No.	Description	Hub Color	
26500	18G x 1"	Pink	
26501	20G x 1"	Yellow	
26502	20G x 1 1/2"	Yellow	
26503	21G x 1"	Deep Green	
26504	21G x 1 1/2"	Deep Green	
26505	22G x 1"	Black	
26506	22G x 1 1/2"	Black	

100/box 1,000/case



MULTI-SAMPLE LUER ADAPTERS

Catalog No.	Description
26537	Luer Adapter

100/bag 1,000/case



MULTI-SAMPLE LUER HOLDERS

Catalog No.	Description
26530	Holder

100/bag 1,000/case (non-sterile)

MULTI-SAMPLE HOLDER W / PRE-ATTACHED LUER ADAPTERS

Catalog No.	Description
26532	Holder w / adapter

50pcs/box 200/case

BN 26532	Holder w / adapter
----------	--------------------

150/bag 1,500/case (bulk, non-sterile)



INJECTION PLUG

Catalog No.	Description
26539	Injection Plug with cap LL

*Not for use with needles larger than 18 gauge.

50/box 200/case



SPINAL NEEDLES

- ✔ This product does not contain Natural Rubber Latex
- ✔ Convenient blister packaging
- ✔ Non-Pyrogenic
- ✔ Transparent plastic hub
- ✔ Siliconized for smooth comfortable penetration



‘QUINCKE’ BEVEL SPINAL NEEDLES

Catalog No.	Description	Hub Color
26960	18G x (1.20mm) x 3½”	Pink 
26964	20G x (0.90mm) x 3½”	Yellow 
26965	22G x (0.70mm) x 2¾”	Black 
26967	22G x (0.70mm) x 3½”	Black 
26970	25G x (0.50mm) x 3½”	Orange 

50 / box 100 / case



SAFELET IV CATHETERS



- ✓ Easy dispenser pack
- ✓ Functional design casing permitting one action opening and no rolling
- ✓ Color-coded casing cap allows for easier identification of catheter size
- ✓ Translucent catheter hub and flashback chamber allows for easy detection of blood flashback at vein insertion
- ✓ Radio-opaque catheter
- ✓ Non-Pyrogenic
- ✓ Precision finish FEP catheter assures stable flow and eliminates catheter tip kink during venipuncture
- ✓ Can be connected to syringe by removing filter cap to expose Luer taper end
- ✓ Use of hydrophobic membrane filter eliminates blood leakage
- ✓ Close and smooth contact between catheter tip and inner needle enables safe and smooth venipuncture
- ✓ Crimped smooth engagement between catheter and hub, without use of adhesive, helps eliminate blood residuals and helps reduce blood coagulation or hemolysis

IV CATHETERS

Stainless Cannula

Item No.	Size	Gauge	Outside Diameter (mm)	Total Length (mm)	Hub Color	
26726	14G x 1 ¼"	14G	1.63	51.0	Orange	■
26729	14G x 2"	14G	1.63	67.0	Orange	■
26730	16G x 1 ¼"	16G	1.2	50.0	Medium Grey	■
26731	16G x 2"	16G	1.2	66.0	Medium Grey	■
26735	18G x 1 ¼"	18G	0.9	47.0	Green	■
26736	18G x 1 ½"	18G	0.9	54.0	Green	■
26738	18G x 2"	18G	0.9	66.0	Green	■
26741	20G x 1"	20G	0.71	41.0	Pink	■
26742	20G x 1 ¼"	20G	0.71	48.0	Pink	■
26743	20G x 2"	20G	0.71	66.0	Pink	■
26746	22G x 1"	22G	0.56	40.0	Blue	■
26751	24G x ¾"	24G	0.41	34.0	Yellow	■

50/box 200/case


SURGICAL BLADES

DISPOSABLE SCALPELS

- ✔ This product does not contain Natural rubber latex
- ✔ Made from premium quality Carbon Steel and Stainless Steel
- ✔ Computer controlled grinding, sharpening and polishing, ensures consistency smooth incision properties
- ✔ Sterilized by GAMMA radiation 2.5 - RADS
- ✔ A Vapor Corrosion Inhibitor (VCI) paper lining, leakage testing, and an aluminum foil peel pack help ensure quality and integrity when ready for use
- ✔ Blade rigidity is tested to eliminate blade breakage
- ✔ Blade are ultrasonically cleaned prior to packing



STAINLESS STEEL SURGICAL BLADES		CARBON STEEL SURGICAL BLADES		STAINLESS STEEL DISPOSABLE SCALPELS	
Catalog No.	Size	Catalog No.	Size	Catalog No.	Size
29500	10	29526	11	29550	10
29502	11	29532	15	29552	11
29508	15	100/box 1,000/case		29556	15
29515	22			29559	22
100/box 1,000/case				10/box 100/case	



IV ADMINISTRATION SETS VENTED/NON-VENTED

- ✓ This product does not contain Natural Rubber Latex
- ✓ Vented/Non Vented
- ✓ Sterile
- ✓ Roller clamp
- ✓ Non-Toxic
- ✓ Occlusion clamp
- ✓ Non-Pyrogen fluid path
- ✓ Pouch Packaging

IV ADMINISTRATION SETS / VENTED / NON VENTED

With Y- Injection Site, Vented / Non vented, with Clamp • 78" Tubing


Priming Volume
23mL/Approx
Inside Diameter
3.0 ± 0.05mm
Outside Diameter
4.1 ± 0.10 / -0.05



Catalog No.Description

29081 15 Drop/mL Luer slip
29083 60 Drop/mL Luer slip


50 pcs/case



Catalog No.Description

29082 15 Drop/mL Luer lock
29085 60 Drop/mL Luer lock

50 pcs/case



Catalog No.Description

29086 15 Drop/mL Lock/Slip
27086 15 Drop/mL Lock/Slip with Needle-Free Injection Site

50 pcs/case

IV ADMINISTRATION SETS / VENTED / NON VENTED 105"

Priming Volume
28mL/Approx
Inside Diameter
3.0 ± 0.05mm
Outside Diameter
4.1 ± 0.10 / -0.05




Catalog No.Description

29088 15 Drop/mL Lock/Slip, Vented/ Non vented, 2 Injection Sites, Roller Clamp, 2 Occlusion Clamps

50 pcs/case

FILTERED IV ADMINISTRATION SET VENTED / NON-VENTED

Priming Volume
57mL/Approx
Inside Diameter
3.0 ± 0.05mm
Outside Diameter
4.2 ± 0.05 / -0.05



Catalog No.Description

26085 20 Drop/mL Luer Lock, Vented/ Non vented, Latex Free Injection Sites, Roller Clamp, Occlusion Clamp
84" Tubing

50 pcs/case

IV EXTENSION SETS



- ✓ These devices do not contain natural rubber latex
- ✓ Sterile
- ✓ Non-Toxic
- ✓ Non-Pyrogenic
- ✓ Vented/Non -Vented
- ✓ Roller Clamp
- ✓ Occlusion Clamp

I.V. Extension Sets

Catalog No.	Description	Qty	Priming Volume	I.D.	O.D.
1712	12" IV Extension Set with Luer Slip Connector, DEHP free	100pcs/cs	1.6mL	2.57	4
1720	20" IV Extension Set with Luer Slip Connector, DEHP free	100pcs/cs	2.5mL	2.57	4
1720WY	20" IV Extension Set with Luer Slip Connector, DEHP free with Injection Site	100pcs/cs	2.5mL	2.57	4
1730	30" IV Extension Set with Luer Slip Connector, DEHP free	100pcs/cs	3.6mL	2.52	4
1730WY	30" IV Extension Set with Luer Slip Connector DEHP free with Injection Site	100pcs/cs	3.6mL	2.52	4
1772	72" Universal IV Administration Set with Luer Lock Connector DEHP free	60pcs/cs	2.8mL	1.4	2.5



Tolerances

Inside Diameter $0 \pm 0.05\text{mm}$

Outside Diameter $0.05 \pm 0.15\text{mm}$



BULK NON-STERILE SYRINGES



- ✓ This product does not contain Natural Rubber Latex
- ✓ Easy to read large bold print
- ✓ Crystal Clear Barrels
- ✓ Positive plunger stop to prevent spillage when drawing
- ✓ Smooth plunger motion

1mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
*BN 26048	Luer Slip	4000
*BN 26049	TB Luer Lock	4000
BN 26050	TB Luer Lock	4000

3mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
*BN 26200	Luer Lock	1500
*BN 26201	Luer Slip	1500

5mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
BN 26230	Luer Lock	2000
BN 26231	Luer Slip	2000

5mL syringes are graduated up to 6mL of total graduations

10mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
BN 26265	Luer Lock	1500
BN 26266	Luer Slip	1500

10mL syringes are graduated up to 12mL of total graduations

30mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
BN 26290	Luer Lock	400
BN 26291	Luer Slip	400
BN 26292	Cath Tip	400

30mL syringes are graduated up to 35mL of total graduations



20mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
BN 26280	Luer Lock	1000
BN 26281	Luer Slip	1000

20mL syringes are graduated up to 25mL of total graduations

50mL BULK NS SYRINGE

Catalog No.	Description	Qty Per Case
BN 26300	Luer Lock	400
BN 26301	Luer Slip	400
BN 26302	Cath Tip	400

50mL syringes are graduated up to 60mL of total graduations

Available with printed barrel or non-printed barrel.** Packaging is subject to change.

*Bulk non-sterile products are sold for non-clinical application or as a component (unfinished product) that will require further processing before clinical use.

Customers who purchase these components to be incorporated into finished medical devices for clinical applications are expected to submit their own premarket submissions (510(k), PMA, etc.) to FDA for their final finished device(s).

PRODUCT SPECIFICATIONS

PRODUCT	QUANTITY WEIGHT			INCHES			MILLIMETER		
	PCS	LBS	KGS	LENGTH	WIDTH	HEIGHT	LENGTH	WIDTH	HEIGHT
1/2-1mL U-40 / Insulin	500	4.82	2.19	17	7	6	432	178	152
1mL Insulin	500	5.36	2.43	17	7	6	432	178	152
0.3mL Insulin	500	4.50	2.04	17	7	6	432	178	152
1mL TB 0 Dead - Space	500	5.38	2.44	17	7	7	432	178	178
1mL TB w/Needle	1000	14.50	6.58	23	16	10	584	406	254
1mL TB Only	1000	11.50	5.22	19	13	10	483	330	254
1mL LL Syringe w/Safety Needle	300	5.34	2.42	16	10	8	406	254	203
1mL LL Syringe w/Safety Needle	2000	22.80	10.34	26	17	17	660	432	432
3mL w/Needle	1000	16.45	7.46	28	16	9	711	406	229
3mL Only	1000	14.30	6.49	28	12	11	711	305	279
3mL Syringe w/ Safety Needle	400	9.00	4.08	16	10	16	406	254	406
5mL w/Needle	800	17.60	7.98	27	16	9	686	406	229
5mL Only	800	15.20	6.89	18	15	12	457	381	305
10mL w/Needle	800	26.40	11.97	30	18	12	762	457	305
10mL Only	800	21.80	9.89	30	14	12	762	356	305
Curved Tip Syringe	500	13.10	5.94	20	12	9	508	305	229
20mL Only	500	26.65	12.09	28	16	15	711	406	381
30mL Only	250	17.65	8.01	31	15	10	787	381	254
50mL Only	150	15.60	7.08	18	17	13	457	432	330
Scalp Vein Set Butterfly	500	7.00	3.18	18	12	10	457	305	254
Blood Collection Set	200	3.30	1.50	12	10	8	305	254	203
IV Catheter	1000	3.30	1.50	11	6	6	279	152	152
Hypodermic Needles	2000	6.55	2.97	17	14	5	432	356	127
Multi-Sample Needle	1000	4.60	2.09	16	7	4	406	178	102
Luer Adapter	1000	5.37	2.44	11	7	7	279	178	178
Luer Lock Holders Pre-Attached	1000	4.65	2.11	12	11	10	305	279	254
Luer Holder	1000	6.30	2.86	15	13	12	381	330	305
Spinal Needle	100	1.70	.77	8	6	6.5	203	152	165
IV Administration Set	50	4.31	1.95	15	9	8	381	229	203
IV Extension Set (case of 100)	100	2.58	1.34	10	8	8	254	203	191
IV Extension Set (case of 60)	60	2.15	.97	12	8	4	292	203	102
Surgical Blades	1000	2.74	1.24	8	7	3	203	178	76
Disposable Scalpels	100	2.48	1.12	10	9	4	254	229	102
Injection Plugs w/Cap	200	1.56	0.71	12	8	3	305	203	76
Aluminum Hub Needles	1000	5.82	2.64	10	9	5	254	229	127
Aluminum Hub Blunt Needles	100	1.00	0.45	8	3	5	203	76	51

Specifications vary slightly between sizes.

Packaging subject to change.

Please contact us for information on products not listed here.

EXELINT INTERNATIONAL, CO.

Corporate Headquarters

2500 Santa Fe Avenue
Redondo Beach
California 90278 USA

Phone: (310) 649-0707

Fax: (310) 649-1178

Email: info@exelmed.com

FOR INFO PLEASE CALL

1-800-940-EXEL(3935)

info@exelint.com
www.exelint.com

